

# Pressure Fit System

Designed to Fit **Any Space**

The Handicare Pressure Fit System can be assembled in minutes, making it a great choice for both temporary and permanent portable solutions. The Pressure Fit track exerts equal pressure between the floor and ceiling to stabilize the unit while ensuring ceilings aren't damaged. It also has a bath adapter that enables safe transfers into standard bathing systems. Available in two, three and four-post designs.



## Specifications & Part Numbers



Part Numbers	Retail Price
2-Post: 341500	\$ 4695.00
3-Post: 341505	\$ 7995.00
4-Post: 341510	\$ 10895.00

Ask your dealer for various packaged systems with trolleys.

**Safe Working Load (SWL)**  
440 lb (200 kg)

**Unit Height**  
Expands from 7' - 9'

Distributed By:  
Covenant Access & Mobility Systems  
3939 S. Charleston Pike - Springfield, OH 45502  
888-510-1933 office / 888-510-1932 fax