

Surveyor		Date		Authority	
Professionals name		Client name			
Address			Address		
Telephone			Telephone		
Email			Email		
Estimate copy to					

Model Normal Distributor Discount x .65 (* =pounds sterling) x 1.3=\$USD\$

T1 Height adjustable
Hydraulic *2200

T2 Static height
Fixed height *1344

T3 Height adjustable
Hydraulic Paediatric *2200

T4 Height adjustable
Electric *3607

Coloured liner Colour*: _____
Other than Grey *260


*Colour chart available from customer services. Standard colour is Grey.

Changes to Standard Trolley: Liner & cot sides, up to 200mm *585
 Length (max 2000mm) *350 Stainless Steel Frame *607
 Width (max 1000mm) *350 Dartex envel. side rails *63
 Other measurements *350 Special Flat Cushion *400


Dimensions

Note: only the lettered dimensions are amendable all other dimensions are provided for information only.


All Models



T1, T3 & T4 Adjustable height models



T2 Fixed height model



Additional Changes:
 Extended Hose *55
 Special Liner *400
 Head Cushion *110
 (non-standard size)

Ref	Feature	Standard Dimensions (mm)	Special Dimensions (mm)
A	External width	720	
B	External length	1930 (T3 = 1600)	
C	Internal length (must be 210mm shorter than B)	1720 (T3 = 1390)	
D	Internal width (must be 180mm narrower than A)	570	
E	Total cot side height	180	

Additional Changes List additional changes with change prices below:

Notes:

- If dimension A is increased beyond 800mm the base will automatically be increased to provide stability.
- Maximum permitted external width (A) is 1000mm.
- This shower trolley cannot be used with a plug to retain water.



Please return completed forms to:
 Chiltern Invadex (UK) Limited, Unit 6C Thorpe Drive, Banbury, Oxon OX16 4UZ
 Email: sales@chilterninvadex.co.uk Fax: 01869 365588